## INSTRUCTIONS FOR USE

### CUFFS & HOSES

for Non-Invasive Measuring of Blood Pressure

 **TEMPERATURE LIMITS**

 **PROTECT FROM MOISTURE**

 **LOT CODE**

 **QUANTITY**

 **EXPIRY**

 **MEDICAL DEVICE**

 **EC ACCORDING TO EU MDR 2017/745**

 **CAUTION**

 **NON-STERILE**

 **MANUFACTURER**

 **LATEX FREE**  
Not made with natural rubber latex

 **READ INSTRUCTIONS FOR USE**

 **REFERENCE**

The classification rules may vary from country to country.

According to the EU MDR 2017/745 or the **Australian D Regulations, NIBP Cuffs and Extensions** are products of **Class I**. According to the rules of **FDA (USA) or Canada, the Cuffs and Extensions for the non-invasive measuring of BLOOD PRESSURE** are **Class II**.

- For each type there are matching tube(s) (extensions) used to connect to the appropriate electromedical device (which are available separately).
- No liquid must penetrate the cuff's bladder or its extension. In the eventuality of this happening the cuff and its extension must be discarded. Any penetration of liquids in the bladder could cause damage to the electromedical device and could shorten the time life of the cuff.
- Before using the cuff and its Extension it is imperative to read attentively the electromedical device instruction's manual. It is also necessary to verify that the two ends of the tubes are compatible and it is also advisable to connect and disconnect several times for testing purposes.
- In the case of the cuff being contaminated by blood or any body fluids, it must be discarded immediately
- It must be remembered that the continual cleaning and disinfecting might reduce the life time of this device. We advise you strongly to use the recommended or equivalent products below and not to exceed the indicated time of application or immersion.
- Nissha Medical Technologies reusable cuffs cannot be sterilised.
- Single patient (single use) Nissha Medical Technologies cuffs can be used for one patient only.
- The class and type of protection against electric shocks are closely linked to those of any electromedical devices (sphygmomanometer) to which the cuff and its Extension may be connected to.
- Always read the instruction manual carefully for the particular device and any accessories used with it before the first use of the device.
- Nissha Medical Technologies may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

ISO 13485:2016

HS.NisshaMedical.com/IFU

## CUFFS & HOSES

### for Non-Invasive Measuring of Blood Pressure

Please carefully read the following instructions:

Failure to observe these precautions for use may lead to undesirable medical consequences.

#### Important note:

This insert is designed to provide guidance for the use and handling of Cuffs and Extensions which are designed for the non invasive measuring of Blood Pressure. No reference is made to a specific medical technique. The manufacturer declines any responsibility for problems resulting from improper use of the product.

## I. IDENTIFICATION / FIELD OF APPLICATION

### IDENTIFICATION

Nissha Medical Technologies cuffs and Extensions for the non-invasive measuring of blood pressure can be used with most monitors (sphygmomanometers) with either one or two pneumatic inputs.

Nissha Medical Technologies cuffs are available either with or without bladders.

They come in two types:

- reusable (can be used on a number of patients after the application of normal hygiene procedures for this product)
- for single-patient use (to be used and/or repositioned uniquely on the same patient – not for use on another patient)

**For each type there is one or two matching tube(s) (extensions) used to connect to the appropriate electromedical device (which are available separately).**

### FIELD OF APPLICATION

Nissha Medical Technologies offers a full range of cuffs which can be used in the majority of cases (adults, children, infants, neonates, premature babies) for:

- occasional measuring of the blood pressure
- medium or long-term monitoring of the blood pressure.

If you wish to place any subsequent orders, please use the code number shown either on the Nissha Medical Technologies cuff or its extension or on its packaging;

or: contact Nissha Medical Technologies or visit the company's website at [HS.NisshaMedical.com/IFU](http://HS.NisshaMedical.com/IFU) for further information about this product.

## II. IDENTIFICATION / FIELD OF APPLICATION

### STORAGE

The Nissha Medical Technologies cuff and extensions must be stored in its original packaging in order to avoid any damage which might reduce its service life and/or performances.

The storage conditions are as follows:

- Ambient temperature: 10 to 50 °C- 50 to 122 °F

### PACKAGING

Nissha Medical Technologies reusable cuffs and extensions are packaged singly.


**III. PERFORMANCES / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL INTEGRITY / PNEUMATICS / ALLERGICITY**

**PERFORMANCES / RELIABILITY**


The performances of Nissha Medical Technologies cuffs are bound up both with those of the medical device (sphygmomanometer) to which they are connected and with correct storage, use and maintenance rules mentioned in this instruction manual.

We recommend that you check that the electromedical device and cuff are working properly before you proceed to take any measurements of the patient. You should also ensure that the cuff is correctly adapted to the patient's limb.

**1 OR 2-INPUT REUSABLE CUFFS WITH BLADDER**

**1 OR 2-INPUT REUSABLE CUFFS WITHOUT BLADDER**

**Nissha Medical Technologies cuffs and extensions** are checked on the basis of a protocol drawn up in accordance with the current standards and regulations applicable to these cuffs. The results of the tests have been confirmed by an accredited Laboratory. **(LNE technical report no. H060007)**

**SAFETY**

**Nissha Medical Technologies cuffs and extensions** fulfill the specifications of the relevant current European and/or international standards: **EN IEC 60601-1/ EN IEC 60601-2-30 / AAMI SP10.**

The safety class and degree of protection against electric shocks (**BF, CF**) are defined by the electromedical device to which the cuffs are connected. Always read the instruction manual for the particular electromedical device before application.

Low frequency leakage currents are measured in accordance with the recommendations of the current standards which are applicable to this product and are below the authorised limits.

Cuffs and extensions not connected to the electromedical device are not protected against the penetration of liquids, so we highly recommend fitting a specially designed seal plug to the connecting tube.

**(Please see the commercial NIPB catalogue)**

**No liquid must penetrate the cuff's bladder. In the eventuality of this happening the cuff must be discarded. Any penetration of liquids in the bladder could cause damage to the electromedical device and could shorten the time life of the cuff.**

**COMPATIBILITY**

On the company's website ([www.NisshaMedical.com](http://www.NisshaMedical.com)), Nissha Medical Technologies offers its customers a downloadable document including information on the appliance's compatibility along with technical details about it.

**MECHANICAL / PNEUMATIC INTEGRITY**

In order to retain the mechanical and pneumatic integrity of Nissha Medical Technologies cuffs, the instructions described in this instruction manual must be properly followed.

When it is not in use, a **cuff** and its extension should be stored in a safe place in order to avoid any damage. Any failure to adhere to the application rules set out in this instruction manual may lead to inaccurate measurements.

**ALLERGICITY**

The materials used to manufacture **Nissha Medical Technologies cuffs and extensions** have been subjected to allergicity tests. These tests have not shown the presence of any products which might cause an intolerable allergic reaction.

(Lab report no. CE DPC 034 PF 0040)

**IV. INSTALLATION / USE / MAINTENANCE / HYGIENE**

**Before using the cuff it is imperative to read attentively the electromedical device instruction's manual. It is also necessary to verify that the two ends of the tubes are compatible and it is also advisable to connect and disconnect several times for testing purposes.**

Please follow these instructions for optimum installation and operation of the cuff:


**INSTALLATION**

On the patient's side:

- Choose a cuff which is suitable for the place/site of measuring.
- Place the inside part of the cuff onto the chosen limb taking care to make sure that the vertical arrow line lines up with the patient's artery. (fig. a)
- Direct the cuff's input tube(s) towards the patient's hand and position the cuff so that the cuff tube's input is at a distance from the brachial (or femoral) artery. (fig. b)
- Fasten the cuff using the Velcro strip making sure that the white arrow placed on the outside part of the cuff is between the two straight white lines which are in the inside of the cuff. (fig. c) - *If it does not, choose a larger or smaller cuff that fits correctly* -
- Connect the cuff to the appropriate connecting tube(s).




Example showing the connection of a connecting tube to the cuff:

On the connecting tube side (extension):      On the cuff side:

**Please note:**




- The cuff must not be placed onto a limb where another device (drip, blood transfusion, etc.) is connected or where the skin is damaged.
- As far as possible avoid fastening the cuff too tightly, as to prevent the stoppage of the venous circulation between each inflation.
- For frequent blood pressure measurings, carry out regular checks on the state of the patient’s tissues around the measuring place.
- Make sure when placing the connecting tube that it will not hinder the patient’s movements.
- Please check that there is no leakage from the cuff’s connector.
- The inflatable part of the cuff must entirely surround the chosen limb for the measurement.

**On the side of the electromedical device:**

- Fasten the connecting tube to the appropriate electromedical device (sphygmomanometer).
- Switch on the device and make sure it is working properly (read the electromedical device’s instruction manual).

Example of connecting tube:

**HOW TO USE:**

**General conditions of use:**

The conditions of use **for cuff and its extension** are as follows:

- Ambient temperature: 10 to 50 °C- 50 to 122 °F

**PREVENTIVE MAINTENANCE:**

Periodically check that the cuff, the Extensions and their connectors are in good condition.

Periodically make sure that there is no loose gap in the parts which connect together which could cause a pneumatic leak.

**HYGIENE:**

**In the case of the cuff being contaminated by blood or any body fluids, it must be discarded immediately**

In order to use a **reusable Nissha Medical Technologies cuff** in departments where sterility levels must be adhered to, **Nissha Medical Technologies** highly recommends using the cleaning or disinfection methods described below:

**Note:** For Vitalcheck brand cuffs, remove the inner bladder from the cuff sleeve during cleaning.

RE-PROCESSING INSTRUCTIONS OF MEDICAL REUSABLE NIBP CUFFS (Acc. to EN ISO 17664:2004)

<b>Manufactured by:</b> Nissha Medical Technologies	
<b>Device(s):</b> reusable cuffs with or without bladders for the non-invasive measuring of blood pressure	
<b>CAUTION</b>	Cuffs not connected to the electromedical device are not protected against any penetration of liquids, so we highly recommend fitting a specially designed seal plug to the connecting tube.
<b>Re-processing limits</b>	Nissha Medical Technologies cuffs are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cuff may shorten its life time
<b>INSTRUCTIONS</b>	
<b>In the place of use:</b>	Remove any stains or dirt with a soft disposable cloth
<b>Isolation and transportation:</b>	Remember always to reprocess the device after use.

<b>Cleaning preparation:</b>	PREPARE a cleaning detergent such as Enzol, Cidezyme or Burnishine. Follow the instructions of the manufacturer for proper solution preparation.
<b>Automatic cleaning:</b>	Is strictly not applicable
<b>Manual cleaning:</b>	<ul style="list-style-type: none"> <li>• APPLY liberal amounts of the cleaning detergent on the cuff and for one and half minutes.</li> <li>Please note: If the cuff has to be submerged or soaked in a liquid, the cuff tubes must be capped so that the liquid does not penetrate the cuff tubes.</li> <li>• USE a soft brush or cleaning cloth to remove stains.</li> <li>• REPEAT as is necessary</li> <li>• Afterwards RINSE thoroughly with water.</li> </ul>
<b>Disinfection:</b>	<ul style="list-style-type: none"> <li>• PREPARE a disinfecting solution such as Cidex, Isopropyl Alcohol, Sporicidin or 10% bleach solution in separate containers. Follow the instructions of the manufacturer for proper preparation of the solution.</li> <li>• APPLY disinfecting solution on the cuff until it is fully saturated for not more than 6 minutes.</li> <li>• Afterwards RINSE thoroughly with water.</li> </ul>
<b>Drying:</b>	Allow cuff to dry normally without a machine before use.
<b>Maintenance, control and tests:</b>	Please check visually that the cuff is in working order
<b>Preservation:</b>	Please read the conditions of storage in this document
<b>Complementary information:</b>	The above methods have been validated and are strongly advisable. Particularly we advise you strongly to use the recommended or equivalent products above and not to exceed the indicated time of application or immersion, otherwise the life time of the cuff might be reduced.
<b>Contact the manufacturer:</b>	See below or visit the web site <a href="http://www.NisshaMedical.com">www.NisshaMedical.com</a>

**Please note:**

**It must be remembered that the continual cleaning and disinfecting might reduce the life time of this device. We advise you strongly to use the recommended or equivalent products below and not to exceed the indicated time of application or immersion.**

**Nissha Medical Technologies reusable cuffs cannot be sterilised**

**ALLERGICITY:**

**Nissha Medical Technologies cuffs and their extensions do not contain latex.**

**V. GUARANTEE / LIABILITY**

Any unused cuffs and extension kept in their original packaging and which have not been subjected to any visible damage are guaranteed by Nissha Medical Technologies for one year.

Nissha Medical Technologies guarantees the device's conformity to the specifications of the current safety and performance standards applicable to it.

**PLEASE NOTE:**

The class and type of protection against electric shocks are closely linked to those of any electromedical devices (sphygmomanometer) to which the cuff may be connected to.

Always read the instruction manual carefully for the particular device and any accessories used with it before the first use of the device.

Nissha Medical Technologies may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

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 Fabricant

 Sans latex  
Not made with natural rubber latex

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