CLARAVUE®

INSTRUCTIONS FOR USE CLARAVUE® ECG CABLE

	TEMPERATURE LIMITS		ROHS COMPLIANT
	KEEP AWAY FROM DIRECT SUNLIGHT		NON-STERILE
	PROTECT FROM MOISTURE		MANUFACTURER
	LOT CODE		LATEX FREE
	QUANTITY		PHT FREE
	EXPIRY		READ INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED		REFERENCE
	MEDICAL DEVICE		EC ACCORDING TO EU MDR 2017/745
	SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT		

-  The classification rules may vary from country to country. According to European Regulation MDR 2017/745 (Annex VIII) or Australian D Regulations, ECG Cables for CLARAVUE® are products of Class I. According to the rules of FDA (USA) or Canada, ECG cables for CLARAVUE® are Class II.
- Please always remember to read the instruction manual for the electromedical device before using the ECG cable for CLARAVUE®.
- Please keep the ECG cable at a distance as far as possible from any sources of electromagnetic radiation.
- Please in the operating room make sure that the ECG cable for CLARAVUE® (wholly or partly) is outside of the operating area in order to reduce the risk of burning the patient while the electrosurgical unit is in operation.
- In the operating rooms only use the specific high frequency protection ECG cable for CLARAVUE®.
- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- The NISSHA MEDICAL TECHNOLOGIES ECG cable for CLARAVUE® Prewired Electrodes has not been designed to be used in MRI environment.
- If you need to take a cable out of the department for maintenance or verification, this becomes the responsibility of the user department to clean and disinfect the product before shipping or transportation.
- Do not allow any liquid to get inside the various connectors on the interface cable.
- The class and type of protection (BF, CF) against electric shocks are defined by the type of electromedical device to which the cable ECG for CLARAVUE® is connected.
- Always read the instruction manual for the particular device and any accessories to be used before placing the device into service.
- NISSHA MEDICAL TECHNOLOGIES may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

HS.NisshaMedical.com/IFU

CLARAVUE® ECG CABLES

Please read carefully the following instructions

Failure to observe these precautions for use may lead to undesirable medical consequences.

Important note:

This insert is designed to provide guidance for the use and handling of NISSHA MEDICAL TECHNOLOGIES ECG Cable for CLARAVUE® ECG Electrodes Strip which is designed for the measurement of the electrocardiogram in order to diagnose or monitor.

No reference is made to a specific medical technique. The manufacturer declines any responsibility for problems resulting from improper use of the product.

I. IDENTIFICATION FIELD OF APPLICATION

IDENTIFICATION:

The NISSHA MEDICAL TECHNOLOGIES ECG Cable for CLARAVUE® ECG electrodes strip is a crucial part of the system used to monitor or diagnose

Harness connection to interface cable & example of interface cable
Harness connection to interface cable & example of interface cable

the cardiac electrical signal (ECG).

It allows a safe connection between

the electrodes on the patient and the electromedical device.

It is made of:

- a trunk cable fitted with two specific connectors:
 - one of them connects the electromedical device
 - and the other contains either one or two plugs on the patient's side for one or two CLARAVUE® strips,

FIELD OF APPLICATION:

NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip can be used on any kind of patient (adult, infant, paediatric and neonate) for both long and medium term ECG monitoring, resting ECG or for stress tests, on an outpatient basis, in the operating rooms, in Emergency sectors or in any other discipline requiring the measuring of an ECG.

It must be highly noted that the NISSHA MEDICAL TECHNOLOGIES ECG Cable for CLARAVUE® ECG electrodes strips does not contain almost any magnetizable material (that is traces of nickel). Their utilisation in MRI room is not recommended.

ANTIBACTERIAL QUALITIES

The envelopes of the medical bioactive NISSHA MEDICAL TECHNOLOGIES cables for CLARAVUE® carry an additional chemical compost of non toxic material which gives them a slightly acid surface like the human skin pH neighbour 5. Almost all the pathogenic agents responsible for the nosocomial sicknesses are destroyed in this environment.

NISSHA MEDICAL TECHNOLOGIES offers a full range of ECG cables for CLARAVUE® in its catalogue which can be viewed and downloaded on the company's website at: nisshamedical.com

If you wish to place any subsequent orders, please use the code number shown either on the ECG cable for CLARAVUE® or on its packaging.

If you require any further information about this product, please contact NISSHA MEDICAL TECHNOLOGIES SAS or visit the company's website at nisshamedical.com

II – STORAGE / PACKAGING

STORAGE:

The storage conditions for NISSHA MEDICAL TECHNOLOGIES ECG cables for CLARAVUE® are as follows:

- Ambient temperature: -40 to +70° C (-40° to +158° F)
- Relative humidity: 10 to 95% (with no condensation)
- Atmospheric pressure: 500 to 1060 hPa (mb) (7.25 to 15.37 psi/14.8 to 31.39 inHg)

PACKAGING:

Depending on the type (one-piece or detachable), **NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip** are packaged either singly or in separate units (trunk cable, leadwires).

When it is not in use, the **ECG cable for CLARAVUE®** must be stored in its original packaging in order to prevent any damage which might reduce its life time, performances and/or safety level.

2.3 Packaging

Each set / harness is packaged in one OPP/PE laminated pouch	Shipping cartons contain 12 boxes
Depending on the pouch size, 10 to 50 are supplied per box	Sets are supplied non-sterile

In order to be able to identify more easily each ECG set reference is associated to a coloured graphic logo.

In case of several different CLARAVUE® types being used in the same department, the coloured logo will help to save time and avoid errors.

CLARAVUE® ECG electrode strips are referred to as “**devices**” in the subsequent chapters of this instruction manual.

III – PERFORMANCES / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL INTEGRITY / SYMBOLS / ALLERGICITY

PERFORMANCES / RELIABILITY:

NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip are inspected both during and at the end of the manufacturing process according to technical protocols drawn up in line with current standards and regulations of this type of product.

The final results of the tests carried out on a group of representative ECG cables for CLARAVUE® have been confirmed by an accredited Laboratory. (LNE technical report no. J021341.)

They have also undergone clinical testing and assessments.

SAFETY:

NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip are designed and manufactured in accordance with the general and special specifications of the relevant current European and/or international standards:

International standards IEC 60601-1 & 60601-2-25 / 60601-2-27

American standard AAMI EC 53

NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip are part of the “applied part” to the patient as defined by the EN IEC 60601-1 International safety standard.

The **ECG cables’ safety class, type of protection (BF, CF) and degree of protection** against electrical shocks are closely linked to those of the electromedical device to which it may be connected to.

- **Please always remember to read the instruction manual for the electromedical device before using the ECG cable for CLARAVUE®.**

- **Please keep the ECG cable for CLARAVUE® at a distance as far as possible from any sources of electromagnetic radiation.**

- **The NISSHA MEDICAL TECHNOLOGIES ECG cable for CLARAVUE® has not been designed to be used in RMI environment.**

- Please in the operating room make sure that the **ECG cable for CLARAVUE®** (wholly or partly) is outside of the operating area

in order to reduce the risk of burning the patient while the electrosurgical unit is in operation.
- In the operating rooms only use the specific high frequency protection ECG cable for CLARAVUE®.

NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip are designed to withstand repeated defibrillation shocks.

PLEASE NOTE:

NISSHA MEDICAL TECHNOLOGIES may not be held liable for any incidents which might occur to the patient or to medical staff caused by the presence of dangerous electrical currents coming from the electromedical device.

COMPATIBILITY:

On the company's website (nisshamedical.com), NISSHA MEDICAL TECHNOLOGIES offers its customers a downloadable document including information on the device's compatibility along with technical details about it.

MECHANICAL and ELECTRICAL INTEGRITY:

The NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip (conductors, connectors, cable, etc.) offer good mechanical resistance to traction and to flexion/bending in order to reduce the risk of damage during use, NISSHA MEDICAL TECHNOLOGIES SAS has used materials offering high quality and excellent reliability.

The moulded or assembled connectors are fitted with soft inserts designed to minimise the risk of the cable breaking at this point.

The connector contacts are treated to reduce contact resistance. (Test report no. CPB1097 – CPB1897/CPB2897)

The service life of an ECG cable depends on a large number of parameters. e.g.

- How many times it has been used
- At all times be sure to follow the safety instructions which you will find in the “PREVENTIVE MAINTENANCE” section below
- Being kept in a good state of hygiene

A regular visual and electrical inspection of the ECG cable will determine whether or not it needs to be replaced. (See preventive maintenance)

SAFETY SYMBOLS:

Explanation of the symbols used on the electromedical device:

Shows that the ECG cables and accessories are fitted with special protection against electric shocks (including admissible leakage currents) and defibrillation shocks.

The class and type of protection (BF, CF) against electrical shocks are defined by the type of electromedical device to which the NISSHA MEDICAL TECHNOLOGIES ECG cable CLARAVUE® is connected to.

ALLERGICITY:

The materials which come into contact with the patient used to manufacture NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip have been subjected to allergicity tests. These tests have not shown the presence of any products which may lead to any intolerable allergic reaction.

IV – INSTALLATION / USE / MAINTENANCE / HYGIENE / STERILISATION

Please follow these instructions for optimum installation and operation of the NISSHA MEDICAL TECHNOLOGIES ECG cable CLARAVUE®:


(See also the instruction manual for the particular electromedical device and accessories)

On the patient side


After having proceeded according to the recommendations described in the specific IFU on the placing of the CLARAVUE® strips/

1. Connect the IP cable plugs to the appropriate CLARAVUE® strips.

PLEASE NOTE:
The colour of the plugs of the electrodes strips corresponds to the colour of the plugs of the receiving cables.

2. Fasten the NISSHA MEDICAL TECHNOLOGIES ECG cable for CLARAVUE® nearest to the patient using the appropriate fastening system located on the lower side of the cable.
3. Connect the plug of the NISSHA MEDICAL TECHNOLOGIES ECG cable fo CLARAVUE® to the appropriate electromedical device. Plug examples:

4. Start the electromedical device (see the instruction manual for the device).

HOW TO USE:

GENERAL CONDITIONS:

The conditions for using the NISSHA MEDICAL TECHNOLOGIES ECG cable for CLARAVUE® are as follows:

- Ambient temperature: +10 to +40°C (+50 to +104°F)
- Relative humidity: 30 to 50 %
- Atmospheric pressure: 700 to 1060 hPa (10 to 15.4 psi or 20.7 to 31.3 inHg)

SPECIAL CONDITIONS:

- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- The NISSHA MEDICAL TECHNOLOGIES ECG cable for CLARAVUE® has not been designed to be used in RMI environment.

PREVENTIVE MAINTENANCE:

- If you need to take a cable out of the department for maintenance or verification, this becomes the responsibility of the user department to clean and disinfect the product before shipping or transportation.
- Make sure that the insulating materials (protective sheath on the connecting cable, interface cable, connectors) have not suffered any mechanical damage.
- Check the electrical continuity of the conductors monthly using an NISSHA MEDICAL TECHNOLOGIES physiological signal simulator.
 (Contact NISSHA MEDICAL TECHNOLOGIES SAS to find out which is the most suitable tester for your needs)

CORRECTIVE MAINTENANCE:

- There is no corrective maintenance for this product. If it is damaged it is to be discarded immediately.

HYGIENE:

Clean the ECG cable and its leadwires regularly using the method described below.

PLEASE NOTE:

Do not allow any liquid to get inside the various connectors on the interface cable.
 Never sterilise the cables in an "autoclave" using boiling water or steam.

When using the NISSHA MEDICAL TECHNOLOGIES ECG cable in departments where a level of sterility must be adhered to, you must proceed to use the cleaning and/or disinfection methods described below.

Manufactured by: NMT SAS	Method
Device(s): ECG cables	
CAUTION	When cleaning or disinfecting the cables make sure that the plugs at the ends of the cables do not become immersed in any liquid to avoid any electrical problem.
Re-processing limits	NISSHA MEDICAL TECHNOLOGIES ECG Cables for CLARAVUE® ECG electrodes strip are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its life time
INSTRUCTIONS	
In the place of use:	Remove any stains or dirt with a soft disposable cloth
Isolation and transportation:	Remember allways to reprocess the device after use.
Cleaning preparation:	No specific requirement.
Automatic cleaning:	Is strictly not applicable
Manual cleaning:	<ul style="list-style-type: none"> • Clean the NISSHA MEDICAL TECHNOLOGIES ECG cable and its connecting cable (including the connector) with a cloth soaked in soapy water
Disinfection:	<p>Method A (Taken from the study by the ANIOS Laboratory, no. 6416.94/0387)</p> <ul style="list-style-type: none"> - fill a dip tank with a 0.5% HEXANIOS G+R solution - partially immerse the cable keeping the ends dry in order to avoid any electrical problems on the connectors - allow a contact time of 15 minutes - rinse the cable - dry the cable with absorbent paper - during the immersion rub the ends of the cable using a LINGET ANIOS wipe. <p>(Replace the solution every 48 hours)</p> <p>Method B (new cold method) (Taken from the study by the ANIOS Laboratory, no. 14496.02/052)</p> <ul style="list-style-type: none"> - fill a dip tank with an ANIOXYDE 1000 activated preparation - but do not use before 30 minutes - partially immerse the cable protecting the ends in order to avoid any electrical problems on the connectors - allow a contact time of 10 to 30 minutes depending upon the required level of disinfection - regularly check the peracetic acid content using a test strip - and afterwards rinse with ANIOS mains water (pH = 7.3; TH = 48°F)
Drying:	Dry the cable with absorbent paper
Maintenance, control and tests:	Please check visually that the ECG cable is in working order
Conditioning:	

Conditioning:	Please follow the standard protocol of packaging of your institution before sterilisation.
Sterilization:	Please follow the standard protocol of packaging of your institution before sterilization.
Please Note:	Never sterilize the cables and accessories in an "autoclave" using boiling water or steam.
Preservation:	Please read the conditions of storage in this document
Complementary information:	The above methods have been validated and are strongly advisable. Particularly we advise you strongly to use the recommended or equivalent products above and not to exceed the indicated time of application or immersion; otherwise the life time of the cuff might be reduced.
Contact the manufacturer:	See below or visit the web site nisshamedical.com

V – GUARANTEE / LIABILITY

Any unused NISSHA MEDICAL TECHNOLOGIES ECG cable kept in its original packaging and which has not suffered any visible damage is guaranteed by NISSHA MEDICAL TECHNOLOGIES SAS for one year against hidden defects.

NISSHA MEDICAL TECHNOLOGIES SAS guarantees that the device complies with the specifications of the current safety and performance standards applicable to the ECG cables.

 PLEASE NOTE:

The class and type of protection (BF, CF) against electric shocks are defined by the type of electromedical device to which the cable ECG is connected.
 Always read the instruction manual for the particular device and any accessories to be used before placing the device into service.
 NISSHA MEDICAL TECHNOLOGIES SAS may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

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