



Securement Products for Respiratory Care

Securing cannulas and tubing is essential to delivering the needed respiratory therapy. Nissha Medical Technologies prioritizes protecting young patients' skin by offering a variety of sensitive securement devices.

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

CPAP CANNULAIDE® SECUREMENT DEVICE

- For nCPAP therapy where a nasal seal is intended
- Works well in humid conditions, even up to 24 hours of continuous use
- Transparent material allows for easier visualization of skin integrity and repositioning if needed
- Hook and loop system prevents the cannula from disconnecting from patient movement
- Made with soft silicone material to avoid skin irritation and align with NICU skin care guidelines

Part #	Ref	Product Description	QTY
700000588	CA100-0-25	Size 0, for Babies less than 700g	25/case
700000589	CA101-0-25	Size 1, for Babies 700g -1250g	25/case
700000590	CA102-0-25	Size 2, for Babies 1251g-2000g	25/case
700000591	CA103-0-25	Size 3, for Babies 2001g-3000g	25/case
700000592	CA104-0-25	Size 4, for Babies Over 3000g	25/case
700000593	CA110-0-1	Variety Pack: 3x CA100, 7x CA101, 10x CA102, 3x CA103, 2xCA104	25/case

***CPAP Cannula Securement Device is not for use with low or high flow oxygen cannulas, nCPAP masks or endotracheal tubes.**



Enhancing the Seal


Pre-cut and available to fit different patient weights and cannula prong sizes.

Less pressure on tender young skin

The size of babies' nares varies over time and as babies age, so the limited, fixed sizes of cannula prongs will often not fit. When prongs don't fit, using a larger size or pushing into the nose is required for a seal. Measuring the force required to create CPAP pressure shows this clearly.

The CPAP Cannulaide Securement® Device enhances a seal without having to press into the baby's nares, with just 1oz of force on average.


Average force into the nares



MINI WHISKERS™ SECUREMENT DEVICE - Premature to Neonatal Patients

- Ideal for premature infants less than 2000g
- Provides gentle adhesion for up to 24 hours
- A gentle securement system for nasal cannulas and feeding tubes
- Velcro strips help keep tubing in position, reduce skin friction, and allow for easy adjustments when needed
- A case contains 20 pouches

Part #	Ref	Product Description	QTY
700000594	MW100-0-120	Exclusively for early-premature babies less than 2000g	120/case



MW100

STICKY WHISKERS™ SECUREMENT DEVICE - Neonatal to Pediatric Patients

- A gentle securement system for nasal cannulas and feeding tubes
- Hydrocolloid material designed to provide gentle adhesion for up to 24 hours
- Hook and loop strips help keep tubing in position, reduce skin friction, and allow for easy repositioning of the cannula

Part #	Ref	Product Description	QTY
700000595	SW102-0-25	One size fits all for infants	25/case



SW102


TENDER GRIPS - Pediatric to Adult Patients

- Soft, skin-toned acrylic material designed to hold cannula tubing in position, generally used around the patient's cheeks
- Internal flap allows cannula adjustment without removing the securement from the patient's skin
- Non-toxic, hypoallergenic, medical grade material
- Moisture-resistant
- Suitable for home, outpatient, extended care, transport and hospital environments
- Provides gentle adhesion for up to 24 hours

Part #	Ref	Product Description	QTY
700000596	1005-0-25	Tender Grips, skin fixation system, adult	25 pairs/case
700000597	1006-0-100	Tender Grips, skin fixation system, adult	Roll of 100
700000598	1007-0-25	Tender Grips, skin fixation system, infant	25/case
700000599	1008-0-25	Tender Grips, skin fixation system, newborn	25/case
700000600	1008-0-10	Tender Grips, skin fixation system, newborn	10/case
700000601	1007-0-10	Tender Grips, skin fixation system, infant	10/case



1005



1007



1008

General Purpose Securement

HYDROCOLLOID STRIPS - All Patients

- Hydrocolloid material that provides gentle adhesion to the skin for 24 hours
- Ideal for securing feeding tube and EKG wires
- Suitable for all patients
- Pre-cut on sheets and ready to apply

Part #	Ref	Product Description	QTY
70000606	1217-0-25	Adhesive Hydrocolloid Strips	25/case



1217

The OEM trademarks identified herein are the trademarks of the respective OEMs, and not of Nissha Medical Technologies. Nissha Medical Technologies disclaims any affiliation, connection, or association between its products and those of the respective OEMs.