

# ClearScan™

Radiolucent Electrodes for ECG Monitoring

## FEATURES



### RADIOLUCENT

Carbon-filled Ag/AgCL sensor is radiolucent for unobstructed imaging.



### HYPOALLERGENIC

Hypoallergenic adhesive and not made with natural rubber latex.



### PACKAGING OPTIONS

Available in economically packaged, no-waste single-use packs or bulk packs.

For use in the Cath Lab, EP, MRI and Radiology or wherever imaging occurs.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

### INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables

MR Conditional



Latex-free



Wet Gel



Solid Gel







# ClearScan™





Radiolucent Electrodes for ECG Monitoring

Keep your imaging free of obstruction by using Vermed ClearScan™ Electrodes. Designed with a carbon-filled plastic snap that will not interfere with your Nuclear Imaging Systems.





## A10004-5 - 1.5" / 38.1mm - TEARDROP FOAM ELECTRODE WITH CARBON SNAP

|  | Product # | Per Strip | Per Pouch | Per Box | Per Case | Wet Gel   | MR Conditional  | Duration of Use   |
|---|-----------|-----------|-----------|---------|----------|---|---|---|
|   | A10004-5  | 5         | 5         | 100     | 600      |  |  |  |
|   | A10004-60 | 5         | 60        | n/a     | 600      |   |   | 24 Hours  |





## A10004-5T - 1.625" / 41.275mm - ROUND TAPE ELECTRODE WITH CARBON SNAP

|  | Product #  | Per Strip | Per Pouch | Per Box | Per Case | Wet Gel   | MR Conditional  | Duration of Use   |
|---|------------|-----------|-----------|---------|----------|---|---|---|
|   | A10004-5T  | 5         | 5         | 100     | 600      |  |  |  |
|   | A10004-60T | 5         | 60        | n/a     | 600      |   |   | 24 Hours  |


## A10006-5RT - 1.875" / 47.625mm - TEARDROP FOAM ELECTRODE WITH CARBON SNAP

|  | Product #    | Per Strip | Per Pouch | Per Box | Per Case | Wet Gel   | MR Conditional  | Duration of Use   |
|---|--------------|-----------|-----------|---------|----------|---|---|---|
|   | A10006-5RT   | 5         | 5         | 100     | 600      |  |  |  |
|   | A1006-1-60RT | 1         | 60        | n/a     | 600      |   |   | 24 Hours  |




## A10025-60 - 1.875" / 47.625mm - TEARTROP FOAM ELECTRODE WITH CARBON SNAP

|  | Product # | Per Strip | Per Pouch | Per Box | Per Case | Wet Gel   | MR Conditional  | Duration of Use   |
|---|-----------|-----------|-----------|---------|----------|---|---|---|
|   | A10025-60 | 5         | 60        | n/a     | 600      |  |  |  |

## A10031-60RT - 1.875" / 47.625mm - TEARDROP FOAM ELECTRODE WITH CARBON SNAP

|  | Product #   | Per Strip | Per Pouch | Per Box | Per Case | Solid Gel   | MR Conditional  | Duration of Use   |
|---|-------------|-----------|-----------|---------|----------|---|---|---|
|   | A10031-60RT | 5         | 60        | n/a     | 600      |  |  |  |

## A10032-5NB - 1.875" / 47.625mm - TEARDROP TAPE ELECTRODE WITH CARBON SNAP

|  | Product #  | Per Strip | Per Pouch | Per Box | Per Case | Solid Gel   | MR Conditional  | Duration of Use   |
|---|------------|-----------|-----------|---------|----------|---|---|---|
|   | A10032-5NB | 5         | 5         | n/a     | 600      |  |  |  |