



EMS Electrodes

12 Lead ECG Kit

FEATURES



V1-V2

Small rectangular electrode with two sensors



V3-V6:

Small rectangular shaped electrode for placement in tight areas




RA, LA, RR, LL:

Large round electrodes for superior adhesion

Our single use pouch contains 10 electrodes perfect for EMS diagnostic stress tests.

PACKAGING



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



CASE

FDA / GS1 Compliant Labels




Built with time savings in mind, our single use pouch contains 10 electrodes perfect for EMS diagnostic stress tests.

A10052-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT
10 ELECTRODES/ KIT | 50 KITS / CASE

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.875" x 1.25" / 47.625mm x 31.75mm	Rectangle	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

WHAT'S IN A KIT?



Perforation allows separation when desired

4 2.25" / 57.15mm
Foam | Wet Gel



4 1.875" x 1.25" / 47.625mm x 31.75mm
Foam | Wet Gel



1 1.75" x 1.625" / 44.45mm x 41.275mm
Perforated | Foam | Wet Gel



Perforation allows separation when desired