



Guardian™

ECG Electrodes for Long-Term Monitoring

FEATURES



HYPOALLERGENIC

Gentle, hypoallergenic adhesives, not made with natural rubber latex.



LOW CHLORIDE WET GEL

Low-Chloride wet gel reduces skin irritation allowing long-term use. Wet gel also hydrates the skin and produces fast, clear traces.



UNIVERSAL SNAPS

High quality Ag/AgCL sensor helps produce clear traces.

Vermed Guardian™ ECG Monitoring Electrodes provide optimal ICU/CCU and other long-term applications.

PACKAGING

RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | HS.NisshaMedical.com

Patient Monitoring & Medical Surgical Consumables

Latex Free



Wet Gel



Guardian™

ECG Electrodes for Long-Term Monitoring

Vermed Guardian™ ECG Monitoring Electrodes provide optimal ICU/CCU and other long-term applications. Guardian electrodes are available in soft cloth, breathable tape, or occlusive foam. Low chloride liquid gel and hypoallergenic adhesives reduce the possibility of irritation associated with long-term use.

A10006-3 - 2" / 50.8mm - ROUND FOAM ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10006-3	3	3	75	600		
A10006-60	3	60	n/a	600		72 Hours

A10006-60T - 2" / 50.8mm - ROUND TAPE ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10006-60T	3	60	n/a	600		
						72 Hours

A10007-3 - 2" / 50.8mm - ROUND CLOTH ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10007-3	3	3	75	600		
						72 Hours