




Defibrillator Compatibility

Match Defibrillator Brands to Nissha Medical Part Numbers

BRAND / MODELS OF DEFIBRILLATOR	COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS	
 <p>Physio-Control LifePak - 15 LifePak - 20 / 20e LifePak - 1000</p>	T100AC-Physio T100LOAC-Physio	C100LOAC-Physio C100AC-Physio *Pediatric-Physio
 <p>Zoll Medical R-Series X-Series</p>	T100AC-Zoll T100LOAC-Zoll T100AC-Zoll-10	C100LOAC-Zoll C100AC-Zoll *Pediatric-Zoll
 <p>Cardiac Science PowerHeart AED G3 Plus PowerHeart AED G3 Pro</p>	T100-CS	

*HeartSync pediatric defibrillation electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Please contact Nissha Medical Technologies at 1.800.669.6905 or HS.Info@NisshaMedical.com for any questions regarding electrode compatibility.