



# VersaTrode™ Wet Gel

Adult General Purpose Electrodes

## FEATURES



### WET GEL ELECTRODES

Wet gel provides quick reading and stabilization.



### VARIETY OF SHAPES AND SIZES

Available in multiple shapes and sizes for better nesting of electrodes in close spaces.



### A10097 VERSATRODE™ LEADLOCKING

Eliminate the need to tape down lead wires or replace displaced electrodes saving time and money. Leadlocking electrodes safely secure leadwires, preventing artifact associated with movement.

Versed VersaTrode™ electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables

Latex-free



Wet Gel



# VersaTrode™ Wet Gel

Adult General Purpose Electrodes

VersaTrode™ disposable electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications. Designed to work in all procedures: stress test, long-term monitoring, short-term monitoring or resting EKG. Reduce inventory by choosing one VersaTrode™ as your house electrode.

## A10005 - 1.5" / 38.1mm - TEARDROP FOAM ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10005-3-7%	3	3	n/a	600		72 Hours
A10005-7%NB	5	5	n/a	600		
A10005-60-7%	5	60	n/a	600		

## A10026T-60 - 1.875" / 47.625mm - TEARDROP TAPE ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10026T-60	5	60	n/a	600		72 Hours

## A10096 - 1.5" x 2.25" / 38.1mm x 57.15mm RECT. FOAM ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10096-5	5	5	n/a	600		72 Hours
A10096-30	5	30	n/a	750		

## A10098 - 1.5" x 2.0" / 38.1MM x 50.8mm - OVAL FOAM ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10098-10B	10	10	n/a	900		72 Hours
A10098-5NB	5	5	n/a	550		

## A10097 - 2.25" x 1.75" / 57.15mm x 44.45mm - LEADLOCKING™ ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10097	1	60	n/a	600		72 Hours
A10097-3	3	3	n/a	600		
A10097-5	5	5	n/a	600		
A10097-5NB	5	5	n/a	600		
A10097-7NB	7	7	n/a	700		